

SILVER



Available in Standard Plastic injection-molded dispenser with receptacle for used blades

Coated with a nitride and copolymer finish to provide smooth wrinkle free sections

They are durable, long lasting and ideal for general paraffin and cryostat use

Offer excellent sections and great blade durability

Made in USA



ORDERING INFORMATION

| Blade | Dimension | Pack size |
|--------------|---------------------|-----------|
| Low Profile | 76mm x 8mm x0.25mm | 50 blades |
| High Profile | 76mm x 14mm x0.43mm | 50 blades |

Manufactured By:



C.L. Sturkey Inc.
824 Cumberland Street
Lebanon, PA 17042
<http://www.sturkey.com>

* Product specification and feature may change without priore information

Exclusively Marketed In India By:

Fugen

Fugen Biomed
11, Dharmabhakti complex, Khushalwadi,
Jamalpure, Navsari, Gujarat 396445
Ph.: +91 909 909 5565