

# INFINITY



## Good as the best and something better than over prices brands

- ☑ Available in Standard Plastic injection-molded dispenser with receptacle for used blades
- ☑ Triple Faceted honed edge with a 35 degree cutting angle.
- ☑ Coated with a special Teflon-like substance to consistently provides a smooth and precise edge
- ☑ Offer excellent sections and great blade durability
- ☑ Made in USA



### ORDERING INFORMATION

Blade	Dimension	Pack size
Low Profile	80mm x 8mm x .25mm	50 blades
High Profile	80mm x 14mm x .30mm	50 blades

Manufactured By:



**C.L. Sturkey Inc.**  
824 Cumberland Street  
Lebanon, PA 17042  
<http://www.sturkey.com>

Exclusively Marketed In India By:

**Fugen**

**Fugen Biomed**

11, Dharmabhakti complex, Khushalvadi,  
Jamalpur, Navsari, Gujarat 396445  
**Ph.:** +91 909 909 5565 | **Email:** [contact@fugenbiomed.com](mailto:contact@fugenbiomed.com)