

CALIBRATORS

CBC CAL PLUS CD CAL NEK CAL TECH CAL



	Vial Volume	Close Vial Stability	Open Vial Stability	Analyzer Group
CBC CAL PLUS	3 mL	45 Days	7 Days	Beckman Coulter and other
CD CAL	3 mL	45 Days	7 Days	Abbott
NEK CAL	3.5 mL	45 Days	5 Days	Systemx
TECH CAL	3.5 mL	45 Days	5 Days	Siemens

* All trademark or register trademark are the property of their respective owners

Availability Schedule

New lot is available EVERY MONTH and expire on 5th of the month after next month

Example

DATE OF AVAILABILITY	EXPIRY
25 th January	5 th March
25 th February	5 th April
25 th March	5 th May

Marketed In India By

fugen
B I O M E D

Fugen Biomed
11, Dharmabhakti Complex,
Khushalwadi, Jamalpore,
Navsari-396445, Gujarat

+91 909 909 5565
contact@fugenbiomed.com
www.fugenbiomed.com