
CBC-5D MR

5 PART DIFF CBC CONTROL









CBC-5D MR 5 PART DIFF CBC CONTROL

Independent third party control

Close Vial Stability Open Vial Stability

WBC diff

Level

Vial Volume

75 Days

14 Days 5 Part

3 (L, N, H) 3 mL

- Internationally accepted third party controls and calibrators
- Totally independent and impartial towards Instrument manufacturer
- Deliveries on standing orders or on regular purchase orders
- Pre-determined and reliable delivery schedule of activity and availability
- Rigorous technical assistance
- Peer group quality control program at affordable price for regular user of control

BC 5000	BC 5100	BC 5140	BC 5300	BC 5200
BC 5600	BC 5120	BC 5150	BC 5380	BC 5500
BC 5800	BC 5130	BC 5180	BC 5390	

^{*} All trademark or register trademark are the property of their respective owners

Availability Schedule			
DATE OF AVAILABILITY	EXPIRY		
25 th February	5 th May		
25 th April	5 th July		
25 th June	5 th September		
25 th August	5 th November		
25 th October	5 th January		
25 th December	5 th March		

ISO 15189 standard requirement for the quality control material

- The new version from December 2012 is applied since the 1st of March 2016
- The essential element recommended in the paragraph 5.6.2.2 concerns the use of independent controls in complement or in replacement of controls
 provided by the instrument supplier
- FUGEN BIOMED supply complete range of Hematology controls for different instruments at affordable Price. We can help you to follow the requirement of the ISO 15189 Standard

Independence

We are completely independent from instrument manufacturer and suppliers who supply their own controls conceived for their own systems called "System Controls"

Reliability

Our controls are made out of fresh Human Blood to be similar to patient samples.

Our products are ready to use, stable and homogeneous, they exist in 3 level of concentration to allow to cover the area of reporting results

Service

Understanding critical requirement of IQC for laboratories we regularly maintain schedule of supply and support for all customer query / problem on priority

Manufactured By:





Biotechne 614 McKinley Place NE Minneapolis, MN 55 413 USA Marketed In India Bv:



Fugen Biomed

11, Dharmabhakti complex, Khushalwadi, Jamalpore, Navsari, Gujarat 396445 *Ph.*: +91 909 909 5565 | *Email*: contact@fugenbiomed.com