

# CBC-X

## 5 PART DIFF CBC CONTROL



**Fugen**  
B I O M E D

# CBC-X

## 5 PART DIFF CBC CONTROL

**Independent  
third party  
control**

**Close Vial  
Stability**

**75  
Days**

**Open Vial  
Stability**

**15  
Days**

**WBC diff**

**5  
Part**

**Level**

**3  
( L, N, H)**

**Vial Volume**

**4.5  
mL**

- Internationally accepted third party controls and calibrators
- Totally independent and impartial towards Instrument manufacturer
- Deliveries on standing orders or on regular purchase orders
- Pre-determined and reliable delivery schedule of activity and availability
- Rigorous technical assistance
- Peer group quality control program at affordable price for regular user of control

### Value Sheet Available for Analysers

SYSMEX	XN Series	XE5000	XT 2000 I	XE 2100
	XT 1800 i	XT 2000 i	XS 1000 I	XE 4000 i
MINDRAY	BC 6600	BC 6800		
MACCURA	JASDAQ	DH-510	DH-520	
		DH-F560	DH-F580	
JIASHAN	JASDAQ	DH-510	DH-520	

\* All trademark or register trademark are the property of their respective owners

### Availability Schedule

DATE OF AVAILABILITY	EXPIRY
25 <sup>th</sup> February	5 <sup>th</sup> May
25 <sup>th</sup> April	5 <sup>th</sup> July
25 <sup>th</sup> June	5 <sup>th</sup> September
25 <sup>th</sup> August	5 <sup>th</sup> November
25 <sup>th</sup> October	5 <sup>th</sup> January
25 <sup>th</sup> December	5 <sup>th</sup> March

### ISO 15189 standard requirement for the quality control material

- The new version from December 2012 is applied since the 1<sup>st</sup> of March 2016
- The essential element recommended in the paragraph 5.6.2.2 concerns the use of independent controls in complement or in replacement of controls provided by the instrument supplier
- FUGEN BIOMED supply complete range of Hematology controls for different instruments at affordable Price. We can help you to follow the requirement of the ISO 15189 Standard

#### Independence

We are completely independent from instrument manufacturer and suppliers who supply their own controls conceived for their own systems called "System Controls"

#### Reliability

Our controls are made out of fresh Human Blood to be similar to patient samples.  
Our products are ready to use, stable and homogeneous, they exist in 3 level of concentration to allow to cover the area of reporting results

#### Service

Understanding critical requirement of IQC for laboratories we regularly maintain schedule of supply and support for all customer query / problem on priority

Manufactured By: \_\_\_\_\_

Marketed In India By: \_\_\_\_\_



**Biotechnie**  
614 McKinley Place NE  
Minneapolis, MN 55 413  
USA



**Fugen Biomed**  
11, Dharmabhakti complex, Khushalwadi,  
Jamalpure, Navsari, Gujarat 396445  
**Ph.:** +91 909 909 5565 | **Email:** contact@fugenbiomed.com