

CBC-DH

5 PART DIFF CBC CONTROL



Fugen
BIOMED

CBC-DH

5 PART DIFF CBC CONTROL

**Independent
third party
control**

**Close Vial
Stability**

**75
Days**

**Open Vial
Stability**

**14
Days**

WBC diff

**5
Part**

Level

**3
(L, N, H)**

Vial Volume

**3
mL**

- Internationally accepted third party controls and calibrators
- Totally independent and impartial towards Instrument manufacturer
- Deliveries on standing orders or on regular purchase orders
- Pre-determined and reliable delivery schedule of activity and availability
- Rigorous technical assistance
- Peer group quality control program at affordable price for regular user of control

Value Sheet Available for Analysers

DYMIND	Dh71, DH73, DH76, DH51, DH53, DH56, D1 CRP, D2 CRP, D3 CRP, D5 CRP, D6 CRP, D7 CRP, DF1 CRP, DF3 CRP, DF5 CRP, DF50, DF51, DF52, DF53, DF55, DF56, DM71, DM72, DM78, DM79, DM71X, DM72X, DM78X, DM79X, UN71, UN73, Un76
GENRUI	KT-65, KT-66, KT-6610 VET Vh50
WHEISMEN	AC-610, AC-630
ZYBIO	Z5, Z50, Z51, Z52
PROKAN	PE-7100 / PE-7100VET, PE-9900, PE-7100CRP

* All trademark or register trademark are the property of their respective owners

Availability Schedule

DATE OF AVAILABILITY	EXPIRY
25 th February	5 th May
25 th April	5 th July
25 th June	5 th September
25 th August	5 th November
25 th October	5 th January
25 th December	5 th March

ISO 15189 standard requirement for the quality control material

- The new version from December 2012 is applied since the 1st of March 2016
- The essential element recommended in the paragraph 5.6.2.2 concerns the use of independent controls in complement or in replacement of controls provided by the instrument supplier
- FUGEN BIOMED supply complete range of Hematology controls for different instruments at affordable Price. We can help you to follow the requirement of the ISO 15189 Standard

Independence

We are completely independent from instrument manufacturer and suppliers who supply their own controls conceived for their own systems called "System Controls"

Reliability

Our controls are made out of fresh Human Blood to be similar to patient samples. Our products are ready to use, stable and homogeneous, they exist in 3 level of concentration to allow to cover the area of reporting results

Service

Understanding critical requirement of IQC for laboratories we regularly maintain schedule of supply and support for all customer query / problem on priority

Manufactured By: _____

Marketed In India By: _____



Biotechnie

614 McKinley Place NE
Minneapolis, MN 55 413
USA



Fugen Biomed

11, Dharmabhakti Complex,
Khushalwadi, Jamalpore,
Navsari-396445, Gujarat

+91 909 909 5565

contact@fugenbiomed.com

www.fugenbiomed.com