

Advia Retic Plus R&D Retic-I Plus R&D 3K Retic XERet



| | Vial Volume | Level | Close Vial Stability | Open Vial Stability | Analyzer Group |
|-----------------------------|-------------|-----------|----------------------|---------------------|-----------------|
| Advia Retic Plus | 4 mL | Tri-level | 75 Days | 14 Days | Siemens |
| R&D Retic-I Plus | 4 mL | Tri-level | 75 Days | 14 Days | Beckman Coulter |
| R&D 3K Retic | 3.5 mL | Bi-level | 75 Days | 16 Days | Abbott |
| XERet | 3 mL | Tri-level | 75 Days | 15 Days | Sysmex |

* All trademark or register trademark are the property of their respective owners

Availability Schedule for Advia Retic Plus, R&D Retic-I Plus, R&D 3K Retic

| DATE OF AVAILABILITY | EXPIRY |
|----------------------------|--------------------------|
| 25 th January | 5 th April |
| 25 th March | 5 th June |
| 25 th May | 5 th August |
| 25 th July | 5 th October |
| 25 th September | 5 th December |
| 25 th November | 5 th February |

Availability Schedule for R&D Retic, R&D Retic-I, XERet

| DATE OF AVAILABILITY | EXPIRY |
|---------------------------|---------------------------|
| 25 th February | 5 th May |
| 25 th April | 5 th July |
| 25 th June | 5 th September |
| 25 th August | 5 th November |
| 25 th October | 5 th January |
| 25 th December | 5 th March |

Marketed In India By

fugen
B I O M E D

Fugen Biomed
11, Dharmabhakti Complex,
Khushalwadi, Jamalpore,
Navsari-396445, Gujarat

+91 909 909 5565
contact@fugenbiomed.com
www.fugenbiomed.com