

**SPECIAL
CONTROLS**

SEDRite Plus

Close Vial Stability	Open Vial Stability	Level	Vial Volume
195 Days	30 Days	2 (L1, L2)	4.5 mL

Value Sheet Available for		
Clinical Data Inc. Excyte 10/M, Excyte 40	Diesse Mini-Ves, Ves-matic 10/Easy, Ves-Maric 20	R&R Mechatronics STARRSED
Westergren (Diluted)	Westergren (Undiluted)	Wintrobe

* All trademark or register trademark are the property of their respective owners

Availability Schedule	
AVAILABILITY	EXPIRY
25 th June	5 th January
25 th December	5 th July

Sickle QC

Close Vial Stability	Open Vial Stability	Level	Vial Volume
195 Days	100 Days	2 (Pos, Neg)	2.5 mL

Value Sheet Available for		
Chembio Sickle-STAT	Columbia Calibre Sickle Cell Reagent	Dade Behring Sickle Sol solubility test
J&J SICKLEDEX	Pacific Hemastasis SickleScreen	Streck Sickle-CHECK

* All trademark or register trademark are the property of their respective owners

Availability Schedule	
AVAILABILITY	EXPIRY
25 th March	5 th October
25 th September	5 th April

R&D Glu/HgB Controls

Close Vial Stability	Open Vial Stability	Level	Vial Volume
105 Days	30 days	3 (L1, L2, L3)	1.5 mL

Value Sheet Available for
Hemocue Hb201+
Hemocue B-Hemoglobin
Hemocue Glucose 201

* All trademark or register trademark are the property of their respective owners

FETALtrol

Close Vial Stability	Open Vial Stability	Level	Vial Volume
105 Days	25 thermal cycle	3 (L1, L2, L3)	2 mL

Value Sheet Available for		
Catlag - BD FACScanto-II	Sure-tech	BD FACScalibur

* All trademark or register trademark are the property of their respective owners

Availability Schedule for R&D Glu/HgB Control and FETALtrol		
AVAILABILITY	EXPIRY	
	FETALtrol	R&D Glu/HgB
5 th January	10 th May	15 th April
5 th April	10 th August	15 th July
5 th July	10 th November	15 th October
5 th October	10 th February	15 th January

BODY FLUID CONTROL

Body Fluid I

Close Vial Stability	Open Vial Stability	Level	Vial Volume
75 Days	30 Days	3 (L1, L2, L3)	3 mL

Availability Schedule	
AVAILABILITY	EXPIRY
25 th February	5 th May
25 th April	5 th July
25 th June	5 th September
25 th August	5 th November
25 th October	5 th January
25 th December	5 th March

Value Sheet Available for

Haemocytometer	Coulter LH-700 series	Mindray BC-6800
Coulter Unicel DxH 800/600	Sysmex XT-4000i	Sysmex XN-series

* All trademark or register trademark are the property of their respective owners

IMMUNO PHENOTYPIC CONTROLS

Status Flow

	Vial Volume	Close Vial Stability	Open Vial Stability	Value Sheet Available for
Status Flow	2.5 mL	45 Days	9 thermal cycle	Coulter ELITE Coulter EPICS Coulter EPICS XL BD FACS calibur BD FACS canto II
Status Flow^{LO}	2.5 mL	45 Days	9 thermal cycle	
Status Flow^{PRO}	1.5 mL	45 Days	9 thermal cycle	



Availability Schedule

New lot is available EVERY MONTH and expire on 5th of the month after next month

Example

DATE OF AVAILABILITY	EXPIRY
25 th January	5 th March
25 th February	5 th April
25 th March	5 th May

Marketed In India By

Fugen
B I O M E D

Fugen Biomed
11, Dharmabhakti Complex,
Khushalwadi, Jamalpore,
Navsari-396445, Gujarat

+91 909 909 5565
contact@fugenbiomed.com
www.fugenbiomed.com

*Product specification and feature may change without prior intimation.

Ver: 1.0 07/01/2019