

# QC CARE IMMUNO

Immunology Assayed Control



**fugen**  
BIOMED

# QC CARE IMMUNO

## Immunology Assayed Control

Independent  
third party  
control

### Format

**Lyophilized**

### Close vial stability

**2  
Years** (2-8 °C)

### Open vial stability

**20 Days** (-20 °C)  
**7 Days** (2-8 °C)

### Level

**3**  
( L1, L2, L3)

### Vial Volume

**2  
mL**

### Analytes

|                                |                                     |                                    |                                     |
|--------------------------------|-------------------------------------|------------------------------------|-------------------------------------|
| 17 OH-Progesterone             | Thyroid Stimulating Hormone (TSH)   | Testosterone, Free                 | Luteinizing Hormone (LH)            |
| Folic Acid                     | Osteocalcin                         | Follicle Stimulating Hormone (FSH) | Homocysteine                        |
| Aldosterone                    | Troponin I                          | Vitamin B12                        | Prostatic Acid Phosphatase (PAP)    |
| AlphaFetoprotein               | Estriol, Free (uE3)                 | Myoglobin                          | Progesterone                        |
| Androstenedione                | Cortisol                            | Valproic Acid                      | Prolactin                           |
| CA 125                         | C-Peptide                           | Carbamazepine                      | Prostate Specific Antigen (PSA)     |
| CA 15-3                        | Dehydroepiandrosterone Sulfate      | Digoxin                            | Sex Hormone Binding Globulin (SHBG) |
| CA 19-9                        | Estradiol                           | Digitoxin                          | Phenytoin                           |
| Carcinoembryonic Antigen (CEA) | Ferritin                            | 2 Erythropoietin                   | C Reactive Protein (CRP)            |
| T3, Total                      | Prostate Specific Antigen Free(PSA) | Human Chorionic Gonadotropin (hCG) | Theophylline                        |
| T4, Total                      | T3, Free                            | Human Growth Hormone (hGH)         | Paracetamol                         |
| Testosterone                   | T4, Free                            | Insulin                            | Immunoglobulin E, total (IgE)       |
| Thyroglobulin                  | Osteocalcin                         |                                    |                                     |

### Value Assigned for

- Abbott Architect Systems
- Abbott AxSYM FIA System
- Beckman Access/Dxl CLIA Systems
- Beckman AU Chemistry System
- Beckman CX/LX/DxC Chemistry System
- Biocode Hyclon Lisa
- Biomedis/Biosystems Chemistry Systems
- BioMerieux Vidas FIA System
- Biotecnica Instruments serie BT
- Brahms Kryptor FIA System
- Phadia ImmunoCAP CLIA Systems
- Roche Modular/Cobas6000/8000 Chemistry Systems
- SNIBE Maglumi CLIA Systems
- Tosoh AIA FIA Systems
- Thermo Kone Chemistry Systems
- DiaSorin Liason CLIA Systems
- DIESSE Chorus EIA Systems
- Fujirebio Lumipulse CLIA System
- Horiba ABX
- IL Ilab Systems
- Mindray BS Chemistry Systems
- Ortho Vitros Chemistry Systems
- Roche ModularE/Elecsys CLIA Systems
- Ortho Vitros CLIA System
- Siemens Advia Centaur Systems
- Siemens Advia Chemistry Systems
- Siemens BN ProSpec
- Siemens Dimension Vista
- Siemens Immulite / Immulite1000 & 2000/ 2000 xpi
- Immunodiffusion
- Electrophoresis
- ISE Analyzers
- Nephelometer
- ELISA
- HPLC

\*\*All trademark or register trademark are the property of their respective owners

### Ordering Information

**3 levels x 5 vials x 2 mL**

Manufactured By:

Marketed In India By:



**CARE SRL**  
Via G. Adamoli  
441-16165  
Genova, Italy



**Fugen Biomed**  
11, Dharmabhakti Complex,  
Khushalwadi, Jamalpore,  
Navsari-396445, Gujarat

+91 909 909 5565  
contact@fugenbiomed.com  
www.fugenbiomed.com