



A Quality Urine Strip

URS-10H

Bilirubin

Blood

Glucose

Ketones

Leucocytes

Nitrite

pH

Protein

Specific Gravity / Density

Urobilinogen



- ✓ A Dip and Read test strip
- ✓ Available in single or multiple test format
- ✓ Semi-Quantitative / Qualitative result interpretation through comparison with color charts or with strip reader
- ✓ High quality color chart for accurate visual reading
- ✓ Stable color development and reliable quality of result due to use of high grade raw materials
- ✓ Air tight strip container with desiccant in bottle to prevent reagent deterioration

Marketed In India By:

Fugen
B I O M E D

Fugen Biomed

11, Dharmabhakti Complex,
Khushalwadi, Jamalpure,
Navsari-396445, Gujarat

Customer Care Contact Details

+91 909 909 5565

contact@fugenbiomed.com

www.fugenbiomed.com