

UA101

Semi Auto Urine Strip Analyzer



FAST



ACCURATE



EASY

UA101

Semi Auto Urine Strip Analyser

UA 101 Urine Analyser is a high precision, intelligent instrument which is developed based on modern optics, electronics, computer science and other advanced technologies for clinical testing of urine. It has advantages of easy and quick operation, accurate result and high reproducibility.

TECHNICAL SPECIFICATION

Analyzer Type	Semi-Automatic
Methodology	Reflectance Principle
Detection	Ultra-High brightness cold light source
Display	LCD – 264 x 64 dot matrix, 37 mm x 130 mm
Data Storage	Storage of 2000 test items, freely query by result number
Strip Incubation Time	1 minute
Communication Interface	Standard RS – 232 interface, can connect with urine sediment analysis system, Hospital Management System
Working Environment	Temperature 15°C ~ 35°C, Humidity ≤85%
Reproducibility	CV ≤ 1%
Testing Speed	120 Sample/hour
Printer	Built – in high speed thermal printer
Power Supply	50 VA
Dimension	240mm (L) x 285mm (W) x 140mm(H)
Weight	≤ 3 kg

Marketed By:



Fugen Biomed

11, Dharmabhakti complex,
Khushalwadi, Jamalpore,
Navsari - 396 445, Gujarat

Contact Us:

+91 909 909 5565

contact@fugenbiomed.com

www.fugenbiomed.com